

(CW/Pulse, Up to 20W SM or 10W PM, Benchtop)



The HPSL High Power Single Mode Laser is a user-friendly benchtop unit that delivers up to 20W single mode and 10W polarization Maintain stable laser output. The output is in CW mode or in a variable-modulated mode with a frequency of up to 10kHz. The output is a collimator with a built-in isolator. It has a front power control knob and USB computer interface. An emission switch adds safety.

Features

- Low Cost
- High Reliability
- High Power
- Single Mode
- USB
- Turn-Key Benchtop

Applications

- Lab
- OEM
- Sensor
- Instrumentation

Specifications

Parameters	Min	Typical	Max	Unit	
Operation Wavelength	1055	1060	1072	nm	
Operation Mode		CW / Modulated			
Output Power	2		22	W	
Bam Quality	1.2	1.3	1.3	M2	
Beam Diameter (with collimator)		4		mm	
Spectral Linewidth		2	2	ns	
Polarization Extinction Ratio	18	26	35	dB	
Output Power Adjust Range	10		100	%	
Modulation Range	0		10	kHz	
Output Power Stability (within 48 hr)		± 2	± 5	%	
Operating Temperature	-5		35	°C	
Storage Temperature	-40		85	°C	
Electrical Power Consumption			150	W	
Power Input	110		120	VAC	
Weight		15		lbs	
Computer Interface	USB				

Legal notices: All product information is believed to be accurate and is subject to change without notice. Information contained herein shall legally bind Agiltron only if it is specifically incorporated into the terms and conditions of a sales agreement. Some specific combinations of options may not be available. The user assumes all risks and liability whatsoever in connection with the use of a product or its application.

Rev 04/15/24

P +1 781-935-1200

E sales@photonwares.com

www.agiltron.com



(CW/Pulse, Up to 20W SM or 10W PM, Benchtop)

Operation Manual

- Plug AC power
- Turn ON The Power Switch
- Adjust The Output Power to Minimum by Turning The Knob All Way Counter Clockwise
- Turn On The Emission Switch
- Increase The Out Put Power by Turning The Knob Clockwise
- To Modulator The Laser, Turn On The Modulation Switch at the Back, Input a 0-5V Modulation Signal Via The BNC Connector
- The Laser Can Also Be Controlled By a Computer via The USB/GUI Interface

Mechanical Dimension

*Product dimensions may change without notice. This is sometimes required for non-standard specifications.







(CW/Pulse, Up to 20W SM or 10W PM, Benchtop)

Typical Spectrum

Ordering Information

				2				11
Prefix	Wavelength	Output Power	Mode ^[1]	Modulation	Output Configuration	Power Supply	Interface	
HPSL-	1060nm=1	20W = 2 10W = 1 Special = 0	Random = 1 PMER18dB = 2 PMER25dB = 3 PMER30dB = 4	CW/Modulation= 2 Special = 0	Isolator/Collimator = 1 Isolator/Fiber = 2 Fiber = 3	120-220V = 1	USB = 1 RS232 = 2	

[1]. PMER- Polarization Maintaining Extinction Ratio

Red is Special Order





(CW/Pulse, Up to 20W SM or 10W PM, Benchtop)

USB Command List

Laser Safety

This product meets the appropriate standard in Title 21 of the Code of Federal Regulations (CFR). FDA/CDRH Class 1M laser product. This device has been classified with the FDA/CDRH under accession number 0220191. All versions of this laser are Class 1M laser products, tested according to IEC 60825-1:2007 / EN 60825-1:2007. An additional warning for Class 1M laser products. For diverging beams, this warning shall state that viewing the laser output with certain optical instruments (for example eye loupes, magnifiers, and microscopes) within a distance of 100 mm may pose an eye hazard. For collimated beams, this warning shall state that viewing the laser output with certain instruments designed for use at a distance (for example telescopes and binoculars) may pose an eye hazard.

Wavelength = $1.3/1.5 \mu m$.

Maximum power = 30 mW.



*Caution - Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

*IEC is a registered trademark of the International Electrotechnical Commission.

P +1 781-935-1200

E sales@photonwares.com

www.agiltron.com



(CW/Pulse, Up to 20W SM or 10W PM, Benchtop)

Questions and Answers

Q:



